

ELASTIKOTE® CUSTOM COLOR PROTOCOL: ALL APPLICABLE SUBSTRATES

ElastiKote 100/1000
Bio-Tech 3000/5000

Custom Color Protocol 3 Coat

Preparation

Perform preparation steps per applicable substrate protocol

**ElastiKote White @ 1.0 gal/sq Minimum
If Substrate Is Not Already White**

**ElastiKote Custom Color @ 1.0 gal/sq Minimum
if Second Coat
Otherwise, Apply 1st Coat of Custom Color @1.5 gal/sq**

**ElastiKote Custom Color @ 1.0 gal/sq Minimum
if Third Coat
Otherwise, Apply 2nd Coat of Custom Color @1.5 gal/sq**

Application

Complete preparation steps per applicable substrate protocol.

If this custom color is dark, apply applicable ElastiKote (ElastiKote 100, ElastiKote1000, ElastiKote Bio-Tech 3000, ElastiKote 5000) White roof coating at a minimum of 1.0 gal/sq. the white coating helps prevent subsequent coatings from creating blisters due to extreme roof and ambient temperatures. Apply 2nd and 3rd coats of custom color @ 1.0 gal/sq minimum.

Otherwise, apply initial and 2nd coat of custom color @ 1.5 gal/sq.